



AMXUP DR Upgrade Kit for AMX4 Portables

A lightweight retrofit system to upgrade your AMX portable to use a DR detector quickly and easily.

The AMXUP DR Upgrade Kit converts your AMX4 to work with any wireless DR panel.

30 minutes or less. In less than 30 minutes your technician can be processing new patient orders.

Quick and easy installation. Your existing AMX portable will be upgraded to the latest wireless DR technology and provides your facility with a revitalized state-of-the-art point of patient care.

The kit consists of a durable storage compartment for the detector made of lightweight, fiberglass material that mounts to the front of the AMX body. The detector can easily slide in-and-out of the storage bin and be used with or without an aftermarket protective cover on the detector (*see picture at right*).

A metal mounting bracket specially designed to fit a 14" touchscreen tablet PC is mounted on the vertical column for the technician to drive the software. Industrial strength adhesive mounts are placed on both sides of the tube column to allow the technician to easily move the tablet holder to the left or right side of the column, depending on the patient's position.

TI-BA recommends PerkinElmer XRpad and Vieworks ViVIX wireless 14x17 DR detectors for the AMXUP.

Tablet PC not included with kit. One year warranty on kit components.

