

InnoVet[®] VERSA DIGITAL RADIOGRAPHIC SYSTEM DR



"THE MOST TRUSTED NAME IN VETERINARY RADIOGRAPHY"

INNOVET MEANS INTEGRITY

InnoVet revolutionized veterinary radiography in 1991 with the first x-ray system designed by veterinarians. Since then, InnoVet has manufactured more than 6,000 veterinary x-ray systems, making it the #1 veterinary x-ray manufacturer in the United States.

Add to that, the industry's best warranty and it's no wonder InnoVet is the most trusted name in veterinary radiography.



LOCAL SUPPORT

Capital equipment is a long-term investment that deserves personal attention from a service professional who understands your practice.

Regional InnoVet dealers are factory-trained to install and service your x-ray systems with immediate on-site and remote support. Our trusted support network means image quality will never slow down your practice.

INDUSTRY BEST WARRANTY

Our 7-year limited parts warranty secures your investment today and guarantees your investment value tomorrow.

InnoVet[®]

2901 W. Lawrence Ave. • Chicago, IL 60625
800-972-9776 • 773-588-2444
Fax 773-588-3424
www.innovet4vets.com

SYSTEM SPECIFICATION

DIGITAL DETECTOR

- High quality image resolution: 3.5 LP/mm
- Large 17 in. x 17 in. image field (43 cm x 43 cm)
- 9 megapixel array
- Time to image: less than 10 seconds

SOFTWARE

- User-friendly acquisition and viewing software
- Unlimited viewing licenses
- Share images via DICOM transfer or convert to common formats (.jpg, .bmp, .dcm, .tiff) and burn them to CD/DVD or e-mail
- Integrates with many Practice Management Softwares

TABLE / TUBESTAND

- Encased x-ray tube reduces radiation leakage by over 50% providing extraordinary safety
- Unique column and covers safeguard staff and patients from exposed electrical wires
- Rounded corners prevent hip-level bruising common to veterinary technicians
- Stationary table top: 56 in. x 24 in.

COLLIMATOR

- Radiation leakage meets Federal standards for human medical systems.

HIGH FREQUENCY GENERATOR

- 20 kHz High Frequency technology
- 30 kW power rating
- kVp range: 40 – 125, 1 kV increments
- mA stations: 50-300
- mAs range: 0.2 to 400
- mAs/kVp technique selection
- Automatic generator "Shut-off" circuit
- Adjustable volume, audible exposure signal
- UL Classified

X-RAY TUBE

- 1.0-2.0 mm Focal Spots, 140,000 Heat Units

COMPUTER

- PC with RAID 1 (two 500GB hard drives), 4 GB RAM
- Windows 7 Operating System with XP mode
- 24 inch, 2.3 MP monitor

SYSTEM OPTIONS

Only pay for the extra features you select:

- Anatomical programming upgrade
- 500 mA generator upgrade (300 mA standard)
- Four-way float table top upgrade 67 in. x 29 in.
- Cesium scintillator upgrade